



Proven Solutions for Cost Containment

PMSI Product Knowledge



Sleep Surfaces and Equipment

Pharmacy

**Medical Services
and Equipment**

Settlement Solutions

Sleep Surfaces and Equipment

Product Knowledge

- Bed frames
- Mattresses
- Orthopedic beds
- Cervical pillows
- Bed wedges
- Bed rails
- Traction frames
- Bedroom accessories

Product Knowledge

Standard Hospital Bed Pad – E0251, E0255

- Weight limit: 450lbs
- Bed head/foot and height adjustable by manual cranks
- Can be purchased or rented with/without mattress and rails
 - With mattress (E0255)
 - Without mattress (E0251)



Product Knowledge

Semi-Electric Hospital Bed – E0260

- Weight limit: 450lbs
- Bed height by adjustable manual cranks
- Head/foot of bed by adjustable electric hand control
- Can be purchased or rented with/without mattress and rails



Product Knowledge

Full Electric Hospital Bed – E0265, E0266


- Weight limit: 350lbs
- Pendant control provides motorized positioning of the upper body and knees, and also adjusts bed frame height
- Bed can be purchased or rented separately or as a package including mattress and rails
 - Without mattress (E0265)
 - With mattress (E0266)



Product Knowledge

Bed Mattresses or Pad for treatment and prevention of pressure sores

- Weight limit: 250 – 300lbs
- Available in king, queen, full and twin sizes
- Available in depths of one inch to four inches and as an overlay or replacement mattress
- Foam – E0272
- Gel – E0196
- Air – E0186
- Dry – E0184
- Flotation/water – E0187 and E0198



Product Knowledge

Low Air Loss (LAL) Mattress – E0277

Low air loss and alternating pressure mattress systems are used for prevention as well as treatment for wounds and pressure ulcers. These support surfaces relieve pressure.

- Can be full mattress replacement systems or overlays
- Several different models available
 - Dependant upon patient's wound or stage of the decubitus ulcer



Product Knowledge

Alternating Pressure with Pump - E0181, E0277

- Weight limit: 250lbs
- Must be used on a hospital bed with rails for patient safety

aka air mattress



Product Knowledge

Air Fluidized Replacement by Hill-Rom – E0194

These support devices use a high flow rate of air to fluidize a fine particulate material, such as silicone ceramic glass beads, to produce a support medium that has characteristics very similar to liquid. This produces lower pressure between patient skin and the support medium by uniformly supporting the patient through air fluidization of glass beads or low-air-loss cushions.



Clinitron Rite-Hite AFT



Clinitron At-Home AFT



Clinitron II AFT

Product Knowledge

Lateral Rotation

Lateral rotation reduces the risk of pulmonary complications



Lateral Rotation



Lateral Rotation and Low Airloss



Product Knowledge

Orthopedic Beds/Mattresses – E1399*

The purchase of an orthopedic bed is one of the most important measures to benefit from a good night's sleep, and avoid back pain and muscular discomfort.

- Mattresses come in various levels of comfort
 - Adjustable firmness
 - Plush
 - Pillow top
 - Memory foam
- Two different models available
 - Flat orthopedic
 - Adjustable
- Normal delivery will take three to four weeks to reach the patient
- * *Verify with vendor if a more appropriate HCPCS is available*



Product Knowledge

Tempur-Pedic® Products – E1399

- Many products available including:
 - Orthopedic Mattresses
 - Adjustable beds
 - Overlays
 - Office chairs
 - Wedges
 - Pillows



Adjustable
bed



Overlays



Pillows

Product Knowledge

Cervical Pillows – E0190

- Many different styles available including:
 - Tempur-Pedic®
 - Water-based
 - Eggcrate
 - Roll



Product Knowledge

Bed Wedges – E0190

Bed wedge systems reduce pressure on the spine, relieve muscle tension, increase circulation and reduce stress on the heart.

- Weight limit: 250lbs
- Many different styles available



Product Knowledge

Bed Rails, Full Length – E0310


- Weight limit: 250 – 300lbs
- Available in telescoping for additional cost
- Must be used with hospital bed
- Other options available for patient's home bed



Product Knowledge

Bed Rails, Half Length – E0305

- Weight limit: 250 – 300lbs
- Must be used with hospital bed
- Available in telescoping for additional cost
- Other options available for patient's home bed



Product Knowledge

Rails for Home Bed – E1399, E0305, E0310*

- Weight limit: 250lbs
- Many different styles available
- * *Verify with vendor if a more appropriate HCPCS is available*



Product Knowledge

Trapeze Bar Attached to Bed – E0910

- Weight limit: 250 – 300lbs
- Must be used with hospital bed
- Other options available for patient's home bed



Product Knowledge

Trapeze Bar Free Standing – E0940

- Weight limit: 250 – 300lbs
- Must be used with hospital bed
- Other options available for patient's home bed



Product Knowledge

Fracture Frame

- Weight limit: 300lbs
- Must be used with hospital bed
- Fracture Frame Free Standing – E0930
- Fracture Frame Attached to Bed – E0920



BLG2660Z



YFG-330



RI G2891K



BLG3661u



Product Knowledge

Over Bed Table – E0274

- Weight limit: 250 – 300lbs
- Available in tilt top or with vanity
- Can be used without hospital bed over patient's home bed or recliner

